# Beads- on- Leads (2773007112)



Part Number: 2773007112

73 BEAD ON LEAD

#### **Explanation of Part Numbers:**

- Digits 1 & 2 = Product Class
- − Digits 3 & 4 = Material Grade
- Last digit 1 = Bulk Packed 2 = Taped and Reeled

#### Ferrite suppression beads are supplied assembled on tinned copper wire for automated circuit board assembly.

- Wires are oxygen free high conductivity copper with 100% matte tin plating over a nickel undercoating. The resistance of the wire is 3.5 mOhm for the 22 AWG and 2.2 mOhm for the 20 AWG wire.

### Recommended Soldering Profile

#### Packaging Options:

- Beads- on- leads can be supplied bulk packed. The last digit of bulk packed parts is a "1". Parts with a "2" as the last digit of the part number are supplied taped and reeled per IEC 60286-1 and EIA RS-296- F standards. Taped and reeled parts are supplied 4500 pieces on a 14" reel. Taping details: Component pitch 5 mm. Inside tape spacing 52.5 mm. Tape width 6 mm.
- Our "Bead- on- Lead Suppression Kit" (part number 0199000028) is available for prototype evaluation.


For any bead- on lead requirement not listed here, feel free to contact our customer service group for availability and pricing.

# Catalog Drawing 3D Model

Weight: 0.6 (g)

Dim	mm	mm tol	nominal inch	inch misc.				
A	3.5	±0.25	0.138	_				
В	62	±1.50	2.44					
С	9.5	±0.30	0.374					
D	0.65		0	22 AWG				

Reel Information							
Tape Width	Pitch	Parts 7"	Parts 13"	Parts 14"			
mm	mm	Reel	Reel	Reel			
6	5	_	_	4500			



## **Chart Legend**


+ Test frequency

Typical Impedance	$(\Omega)$
1 MHz	42
5 MHz	93
10 MHz <sup>+</sup>	130
25 MHz <sup>+</sup>	132

Beads- on- leads are controlled for impedances only. Minimum impedance values are specified for the + marked frequencies. The minimum impedance is listed on our catalog drawing.

#### **Catalog Drawing**

The impedance of the 73 & 43 beads- on- leads are measure on the E4990A Impedance Analyzer. The 61 beads- on- leads are tested for impedance on the E4991A / HP4291B Impedance Analyzer.



### **CSV** Download